

The Leader in Falls Prevention


Affordable Monitors To Help Reduce Falls
fallsXTM
08000 232 942 | www.fallsx.co.uk



Bed/Chair Monitor



InfraRed Monitor



TreadNoughtTM

The TreadNoughtTM Floor Sensor is built to last with a rugged construction. See back for details.



NEW!

Sentinel Monitor



£64.95 + VAT
81850

- Sensor Pad or Pull String Magnet Triggers Alarm
- One-Year Bed and Chair Sensor Pads Available
- Connects to Most Nurse Call Systems
- Built-in Cushion Protector Extends Monitor Life
- Powered by 9-volt Battery for Complete Portability
- Optional Mains Power Adapter
- Low Battery Alert
- Sturdy Mounting Clip for Easy Transfer
- Optional Mounting Bracket Available
- Auto-On Feature Eliminates On/Off Switch
- Latex Free
- One-Year Limited Warranty

Pads



1-Year Chair Pad
£29.95 + VAT
83700

1-Year Bed Pad
£42.95 + VAT
82790



Order Online Today at
www.fallsx.co.uk or
Call 0800 232 942

FallsX™ Package - £124.95 + VAT **FX-PACKAGE**

Includes 81850, 82790, 83700 and 0707-578-FX

Our FallsX™ Sentinel package includes the FallsX™ Sentinel monitor, one-year bed pad, one-year chair pad and the universal bed bracket. Save £20.85 over buying each product individually.



InfraRed Monitor

- Dual purpose monitor can be used both to monitor either when patient leaves patient bed or when passes through a doorway
- Beep or Chime alarm tones
- Can be connected to most nurse call systems with optional cable
- Silent alarm option when connected to nurse call system
- Powered by standard 9-volt battery or optional mains power adaptor
- Mounts easily to bed frame, headboard or footboard
- Low-battery indicator
- One-year warranty

Basic Components of the 91170-FX:

- InfraRed Monitor with attached plastic bracket and Velcro strips

Optional Components:

- AC Power Adapter allows for mains power option
- Metal Extension Mounting Bracket, Product No. 0120-125-FX
- Nurse Call Connection Kit, Product No. NCKIT-FX



£79.99 + VAT
91170-FX



Multiple Uses

At FallsX™ we know that your facility needs flexibility, and the FallsX™ InfraRed Monitor can provide just that. Our monitor can be placed at the edge of the bed to make you aware of when someone gets up, or it can be placed in the doorway to let you know when someone exits the room.



TreadNought™

Durable Floor Mats

At TreadNought™ we know that your fall prevention efforts are too important to depend on floor sensor mats that don't last. If you are tired of replacing expensive products that don't last then our TreadNought™ Floor Sensor Mat is the perfect floor mat for you.

The TreadNought™ Floor Sensor Mat is built to last with a rugged construction that far out lasts the competition without added expense.

Our floor sensor mat connects easily to most nurse call systems or remote sounders and can be placed at the bedside, doorway or other locations to activate the nurse call system when a fall risk resident or patient steps onto the sensor.

Most nurse call connections available for immediate shipment for your convenience. Optional anti-slip mesh reduces the potential for slippage on hard surface floors.



Distributed By:

fallsX™

FallsX Freephone: 08000 232 942
www.fallsx.co.uk

REDUCE THE RISK OF FALLS